

CLOSURE PLUS PROTECTION



Clo-Sur^{PLUS} P.A.D.TM

Pressure Applied Dressing

Clo-Sur^{PLUS} P.A.D.TM

Following catheterization, achieving hemostasis with manual compression can be time-consuming and costly. Patients are also at increased risk for catheter-related blood stream infections (CRBSI)¹. The Clo-Sur^{PLUS} P.A.D. antimicrobial dressing accelerates hemostasis and protects against nosocomial infection.



Provides 6 day antimicrobial barrier protection against CRBSI




CLOSURE PLUS PROTECTION

CLOSURE

Non-woven topical pad made with a proprietary polyprolate, a positively-charged material that reacts with negatively-charged blood cells to accelerate the clotting cascade.

Faster hemostasis leads to decreased hold time and time to ambulation.

Can also be used for patients on anticoagulation therapy.



¹ <http://www.safecarecampaign.org/CRBSI.html>

² *Eur Radiol*, 2007 Mar 17(3):693-700.

³ WuXiAppTec Report 867961 May 31, 2011.

PLUS

Unlike other vascular closure devices, the Clo-Sur PLUS P.A.D. has no contraindications and can be used by patients with known sensitivity to bovine-derived materials.

Water soluble for painless and easy removal, which is more comfortable for the patient.

PROTECTION

Provides an antimicrobial barrier for up to 6 days.

Can be used as antimicrobial wrap around indwelling lines to protect from catheter-related blood stream infection (CRBSI).

MICROBIAL ASSAYS TEST RESULTS³

| ORGANISM | 7 DAY INCUBATION | | |
|-----------------|------------------|-------------|-------------|
| | REPLICATE 1 | REPLICATE 2 | REPLICATE 3 |
| A. brasiliensis | No Growth | No Growth | No Growth |
| B. cereus | No Growth | No Growth | No Growth |
| C. albicans | No Growth | No Growth | No Growth |
| E. coli | No Growth | No Growth | No Growth |
| E. faecium | No Growth | No Growth | No Growth |
| P. aeruginosa | No Growth | No Growth | No Growth |
| S. aureus | No Growth | No Growth | No Growth |
| S. pyogenes | No Growth | No Growth | No Growth |

Distributed by:



Order: Catalog Number - **SCV22P-AB** - 4X4cm - packaged 10 per box

Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT www.merit.com

Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 40

Free phone for specific country:

Austria 0800 295374 • Belgium 0800 729 06 (Dutch) 0800 731 72 (French) • Denmark 80 88 00 24 • France 0800 91 60 30 • Germany 0800 182 08 71

Republic of Ireland 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 8184 • Sweden 020 79 24 45 • UK 0800 973 115 • Finland 00800 6374 8633

Manufactured by: Scion Cardio-Vascular, Miami, FL 33186

402389001/B ID 030712