

SPECIFICATIONS


Headway®

Microcatheter with Hydrophilic Coating

1 per box. Includes shaping mandrel and introducer sheath.

Product Name	Description	Product Code	ID (inches)	Tip Shape	Usable Length (cm)	OD Prox./Distal (French)	Tip Markers
Headway17 Advanced 'Soft'	Straight	MC172150S	0.017	STR	150	2.4 / 1.7	2
Headway17 Advanced	Straight	MC172150SX	0.017	STR	150	2.4 / 1.7	2
Headway17 Advanced	Pre-shaped tip 45	MC17215045X	0.017	45	150	2.4 / 1.7	2
Headway17 Advanced	Pre-shaped tip 90	MC17215090X	0.017	90	150	2.4 / 1.7	2
Headway17 Advanced	Pre-shaped tip J	MC172150JX	0.017	J	150	2.4 / 1.7	2
Headway21	Straight	MC212150S	0.021	STR	150	2.5 / 2.0	2
Headway27	Straight	MC272150S	0.027	STR	150	3.1 / 2.6	2
Headway27	Straight	MC272156S	0.027	STR	156	3.1 / 2.6	2

Dead space of Headway® 17 microcatheter is 0.41cc and 0.26cc with MicroVention hub adaptor.



Includes Shaping Mandrel and Introducer Sheath



www.microvention.com

MicroVention, Inc.
 Worldwide Headquarters
 1311 Valencia Avenue
 Tustin, CA 92780 USA
 PH +1.714.247.8000
 PH 1.800.990.8368

MicroVention UK Limited
 Suite 3, The Barracks Building
 10 Cliffords Fort, North Shields
 Tyne and Wear, NE30 1JE UK
 PH +44 (0) 191 258 6777
 F +44 (0) 191 258 5999

MicroVention Europe
 30 bis, rue du Vieil Abrevoir
 78100 Saint-Germain-en-Laye
 France
 PH +33 (1) 39 21 77 46
 F +33 (1) 39 21 16 01

MicroVention Deutschland GmbH
 Hildebrandstr. 4 F
 D-40215, Düsseldorf
 Germany
 PH +49 211 210 798-0
 F +49 211 210 798-29

MICROVENTION, Headway, LVIS, FRED and ERIC are registered trademarks of MicroVention, Inc. PHIL is a trademark of MicroVention, Inc. • Scientific and clinical data related to this document are on file at MicroVention, Inc. Refer to Instructions for Use, contraindications and warnings for additional information. © 2014 MicroVention, Inc.MM304(I) Rev.B 11/14