


VACLOK®

Negative Pressure Syringe



EXPANDING THE POSSIBILITIES



MERTMEDICA
ENDOTEK™

VACLOK®

Negative Pressure Syringe

Merit Endotek's VacLok Negative Pressure Syringe is the perfect tool for creating negative vacuum pressure. With the VacLok, Merit Endotek continues to support clinicians by turning negative pressure into a positive experience.

SYRINGE FEATURES:

HOLDS NEGATIVE PRESSURE

Plunger design holds negative pressure with ease, and won't inadvertently disengage

ERGONOMIC DESIGN

Large wings and skid resistant thumb pad make syringe easy to grip and control

VARIOUS SIZES AVAILABLE

Available in 20ml, 30ml, 60ml

DUAL USE DEVICE


Can be used as negative pressure syringe or as standard syringe

NEGATIVE PRESSURE SETTINGS


Four or six different negative pressure settings (dependent on syringe size)



Step 1: Pull



Step 2: Rotate & Release



VacLok Negative Pressure Syringe

Catalog Number	Size	Pres. Settings	Package
ENDO-VAC20	20ml	4 position	20/box
ENDO-VAC30	30ml	6 position	20/box
ENDO-VAC60	60ml	6 position	20/box



ALIMAXX-ES®
Fully Covered
Esophageal Stent



MAXXWIRE®
Guide Wire



Brighton® Bipolar
Coagulation Probe



BIG60™ Inflation Device



AERO® Fully Covered
Tracheobronchial Stent

