

CORONARY ORBITAL ATHERECTOMY SYSTEM



A dual mechanism of action combines differential sanding and pulsatile forces to safely, effectively and efficiently treat severely calcified






Healthy tissue flexes away from the crown

CROWN SIZING



SPECIFICATIONS

VIPERWIRE ADVANCE WITH FLEX TIP



Tip: Platinum/Tungsten Core: Nitinol with Stainless Steel Support Coil Tip Load at 10 mm: 1.0 g

DIAMONDBACK 360 ORBITAL ATHERECTOMY DEVICE

MODEL NUMBER	CROWN SIZE	LENGTH	QΤΥ
C-2DB-CL125-135	1.25 mm Classic	135 cm	1

VIPERWIRE ADVANCE CORONARY GUIDE WIRE

MODEL NUMBER	DESCRIPTION	LENGTH	γту
GWC-12325LG-FT	0.012"/0.014" Flex Tip	325 cm	5

VIPERSLIDE LUBRICANT

MODEL NUMBER	DESCRIPTION	γту
VPR-SLD2	100 mL Bag	10

OAS PUMP

MODEL NUMBER	DESCRIPTION	QΤΥ
SIP-3000	OAS Pump	1

^{1.} Chambers J, et al. JACC Cardiovasc Interv. 2014;7(5):510-518.
2. Shlofmitz E, et al. Expert Rev Med Devices. 2017;14(11):867-879.

^{3.} CSI data on file.